



**Peninsula Laboratories International, Inc.**

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

**Cat. No. S-3098**

**Thr-Ala-Pro-Arg-Atrial Natriurectic Factor (1-28)  
(Urodilatin) (human, bovine, porcine)**

**Specifications**

Number of Slides: 150  
Sample: Paraffin-embedded or frozen sections  
Shelf Life: One year from production date  
Storage: 4°C  
Assay Time: 2-3 hours following primary incubation

**Kit Components**

- Primary antibody
- FITC-labeled secondary antibody
- Normal rabbit serum
- 5X phosphate buffered saline, pH 7.4 with 0.1% sodium azide
- 3% triton X-100/water
- Glycerin-PBS (9:1)
- Normal goat serum