



Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

Cat. No. S-3161

C-Type Natriuretic Peptide (1-29) (Mouse, Porcine, Rat)

Specifications

Number of Slides: 150
Sample: Paraffin-embedded or frozen sections
Shelf Life: One year from production date
Storage: 4°C
Assay Time: 2-3 hours following primary incubation

Kit Components

- Primary antibody
- FITC-labeled secondary antibody
- Normal rabbit serum
- 5X phosphate buffered saline, pH 7.4 with 0.1% sodium azide
- 3% triton X-100/water
- Glycerin-PBS (9:1)
- Normal goat serum