



## **Peninsula Laboratories International, Inc.**

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

**Cat. No. S-3173**

**rec EGF (human)**

### **Specifications**

Number of Slides: 150  
Sample: Paraffin-embedded or frozen sections  
Shelf Life: One year from production date  
Storage: 4°C  
Assay Time: 2-3 hours following primary incubation

### **Kit Components**

- Primary antibody
- FITC-labeled secondary antibody
- Normal rabbit serum
- 5X phosphate buffered saline, pH 7.4 with 0.1% sodium azide
- 3% triton X-100/water
- Glycerin-PBS (9:1)
- Normal goat serum