



## Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

### **Cat. No. S-4014**

## **C-type Natriuretic Peptide (32-53) (human, porcine, rat)**

### **Specifications**

Number of Slides:	150
Chromogens:	DAB (brown) and AEC (red)
Sample:	Paraffin-embedded or frozen sections
Shelf Life:	One year from production date
Storage:	4°C
Assay Time:	2-3 hours following primary incubation

---

### **Kit Components**

- 30% Hydrogen peroxide Solution
- Blocking Solution
- Primary Antibody
- Biotinylated Secondary Antibody
- SA-HRP Conjugate Solution
- DAB Concentrated Solution
- Substrate Buffer for DAB (ready to use)
- AEC Concentrated Solution
- Concentrated Substrate Buffer for AEC
- Hydrogen Peroxide Concentrate
- Hematoxylin Counterstain
- Mounting Solution