



Peninsula Laboratories International, Inc.

MATERIAL SAFETY DATA SHEET SECTION I: IDENTIFICATION

CHEMICAL NAME C-Type Natriuretic Peptide (32-53) (human, porcine, rat) - Immunohistochemistry Staining Kit, Host: Rabbit
PRODUCT CODE S-4014.0001
CAS NO N/A

SECTION II: INGREDIENTS

PRINCIPLE HAZARDOUS COMPONENTS	PERCENT	THRESHOLD LIMIT VALUE
1% BSA in .1M Sodium phosphate buffer	0.1M	Data Not Available
Thimerosal in Normal Goat serum blocking soln. and Biotinylated AB soln. and HRP-SA conj. soln.	0.2%	Data Not Available
Hydrogen peroxide solution	30%	Data Not Available

SECTION III: PHYSICAL DATA

BOILING POINT N/A	MELTING POINT N/A
VAPOR PRESSURE N/A	SPECIFIC GRAVITY N/A
VAPOR DENSITY (AIR=1) N/A	% VOLATILE BY VOLUME N/A
WATER SOLUBILITY N/A	EVAPORATION RATE (____=1) N/A
APPEARANCE White powder	

SECTION IV: FIRE AND EXPLOSION HAZARD DATA

FLASH POINT N/A
FLAMMABLE LIMITS N/A
EXTINGUISHING MEDIA Appropriate for surrounding fire.
SPECIAL FIRE FIGHTING PROCEDURES Wear protective clothing and self-contained breathing apparatus.
UNUSUAL FIRE & EXPLOSION HAZARDS Burning may produce toxic fumes.

7/2/2014

Peninsula Laboratories International, Inc. 305 Old County Road, San Carlos, CA 94070
Tel: (650) 801-6090 Fax: (650) 595-4071



Peninsula Laboratories International, Inc.

SECTION V: HEALTH HAZARD

THRESHOLD LIMIT VALUE	Data Not Available
HEALTH HAZARDS	May be harmful by inhalation, ingestion, or skin absorption. May cause skin irritation. May cause eye irritation.
FIRST AID PROCEDURES	In case of contact, flush eyes with water. In case of contact, wash skin with soap and water.

SECTION VI: REACTIVITY DATA

STABILITY	UNSTABLE STABLE	X	CONDITIONS TO AVOID Not specific.
INCOMPATIBILITY (MATERIALS TO AVOID) N/A			
HAZARDOUS POLYMERIZATION	MAY OCCUR WILL OCCUR	X	CONDITIONS TO AVOID Not specific
HAZARDOUS DECOMPOSITION PRODUCTS Data Not Available			

SECTION VII: SPILL AND DISPOSAL PROCEDURES

SPILLS
Wear respirator, chemical safety goggle, rubber boots, rubber gloves. Absorb in sand and place in closed container for disposal. Wash spill site.
DISPOSAL
Dispose in accordance with local regulation. Generally, dissolve or mix the material with a combustible solvent and burn in a chemical incinerator equipped with an afterburner and scrubber.

SECTION VIII: PROTECTION INFORMATION

EYE PROTECTION Chemical safety goggles	SKIN PROTECTION Full length lab coat and rubber or vinyl gloves
RESPIRATORY PROTECTION Recommended	VENTILATION RECOMMENDED Mechanical exhaust
OTHER PROTECTION Not Specific	

SECTION IX: SPECIAL PRECAUTIONS

STORAGE AND HANDLING Store in a cool, dry place.
--

The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. It is the users' responsibility to determine the suitability of this information for the adoption of safety precautions as may be necessary. Peninsula Laboratories International, Inc. shall not be held liable for any damage resulting from the handling of the above product.

7/2/2014

Peninsula Laboratories International, Inc. 305 Old County Road, San Carlos, CA 94070
Tel: (650) 801-6090 Fax: (650) 595-4071