



Peninsula Laboratories International, Inc.

MATERIAL SAFETY DATA SHEET

SECTION I: IDENTIFICATION

CHEMICAL NAME Ferritin (human) - Clone mAb 44
PRODUCT CODE T-1320.0200
CLONE mAb 44

SECTION II: INGREDIENTS

COMPONENTS	PERCENT	THRESHOLD LIMIT VALUE
Phosphate buffer saline	pH 7.2	Data Not Available
Bovine serum albumin	10 mg/ml	Data Not Available
Preservative: thimerosal	0.01%	Data Not Available

SECTION III: PHYSICAL DATA

FORMAT	STORAGE
Affinity purified, lyophilized	Original vial: 1 year at 4° - 8°C Stock solution or aliquots thereof: 1 year at -20°C. Avoid repeated thawing and freezing.

For in vitro research only. Caution: this product contains thimerosal, a poisonous and hazardous substance.

7/1/2014

Peninsula Laboratories International, Inc. 305 Old County Road, San Carlos, CA 94070
Tel: (650) 801-6090 Fax: (650) 595-4071