



Peninsula Laboratories International, Inc.

## MATERIAL SAFETY DATA SHEET

### SECTION I: IDENTIFICATION

<b>CHEMICAL NAME</b> F 4/80 Antigen (mouse) - Clone BM8 , FITC labelled
<b>PRODUCT CODE</b> T-2027
<b>CLONE</b> BM8

### SECTION II: INGREDIENTS

COMPONENTS	PERCENT	THRESHOLD LIMIT VALUE
Phosphate buffer saline	pH 7.2	Data Not Available
Bovine serum albumin (BSA)	5 mg/ml	Data Not Available
Preservative: sodium azide	0.09%	Data Not Available

### SECTION III: PHYSICAL DATA

FORMAT	STORAGE
Cell culture supernatant, lyophilized	Original vial: 6 months at 4° - 8°C

For laboratory use and research purposes only. Caution: this product contains sodium azide, a poisonous and hazardous substance.

7/1/2014

Peninsula Laboratories International, Inc. 305 Old County Road, San Carlos, CA 94070  
Tel: (650) 801-6090 Fax: (650) 595-4071