



## Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

### **T-4018** **Rabbit anti-Bradykinin IgG**

**Lot: A07372**

Each vial contains 400ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200µl of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8°C. It is recommended that the reconstituted material be aliquotted and stored at -20°C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

#### **Cross Reaction:**

<b>Peptide:</b>	<b>%:</b>
Bradykinin	100
[Lys]-Bradykinin	100
Biotinyl-Bradykinin	100
[Tyr <sup>01</sup> ]-Bradykinin	100
(Des-Arg <sup>1</sup> )-Bradykinin	100
Bradykinin (4-9)	100
Bradykinin (5-9)	40
Bradykinin (7-9)	0
(Des-Arg <sup>9</sup> )-Bradykinin	0
Bradykinin Potentiator B	0
Angiotensin II	0
Met-Enkephalin	0
Substance P	0
Endothelin-1(human, bovine, dog, mouse, porcine, rat)	0
TRH	0