



Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

T-4036

Rabbit anti-CRF (human, mouse, rat) - IgG

Lot: A07989

Each vial contains 400ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200µl of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8°C. It is recommended that the reconstituted material be aliquotted and stored at -20°C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

Cross Reaction:

Peptide:	%
CRF (human, rat)	100
CRF (porcine)	10
CRF (ovine)	<0.5
CRF (bovine)	<0.5
Prepro CRF (125-151) (human)	0
Sauvagine	5
PACAP-38 (human, mouse, ovine, porcine, rat)	0
LH-RH	N/A
ACTH (1-39) (human)	0
ACTH (7-38) (human)	0
[Arg ⁸]-Vasopressin	N/A
BNP-45 (rat)	0