



Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

T-4079

Rabbit anti-Orexin B (mouse, rat) IgG

Lot: 000509

Each vial contains 400ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200 μ l of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8°C. It is recommended that the reconstituted material be aliquotted and stored at -20°C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

Cross Reaction:

Peptide:	%:
Orexin-B (mouse, rat)	100
Orexin B (human)	100
Orexin-A (human, bovine, mouse, rat)	0
NPY (human, rat)	0
α -MSH	0
Leptin (human)	0
Substance P	0
CRF (human, rat)	0