



## Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

### T-4110

## Rabbit anti-Gastric Inhibitory Peptide (GIP) (human) IgG

**Lot: A08864**

Each vial contains 400ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200µl of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% .

Lyophilized products are stable for long periods when stored at 2-8°C. It is recommended that the reconstituted material be aliquotted and stored at -20°C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

#### Cross Reaction:

Peptide:	%:
GIP (human)	100
GIP (porcine)	100
Gastrin I (human)	0
Big Gastrin I (human)	0
Bombesin	0
Substance P	0
Insulin (human)	0
Amylin (human)	0