



Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

T-4173

Rabbit anti-Prepro-Atrial Natriuretic Factor (ANF) (26-55) (human) IgG

Lot: A01668

Each vial contains 400ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200 μ l of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8°C. It is recommended that the reconstituted material be aliquotted and stored at -20°C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

Cross Reaction:

Peptide:	%:
Prepro-ANF (26-55) (human)	100
γ -ANF (1-25) (human)	100
ANF (1-28) (human)	0
ANF (rat)	0
α -ANF (Chicken)	0
β -ANP (1-28) (dimer, antiparallel) (human)	0
Prepro-ANF (56-92) (human)	0
BNP-32 (human)	0
CNP (32-53) (porcine, human, rat)	0
Atriopeptin III (rat)	0
Insulin (human)	0
rec IGF-1 (human)	0
Cardiodilatin (1-16) (human)	0