



Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

T-4218

Rabbit anti-C-Type Natriuretic Peptide (1-29) (mouse, porcine, rat) - IgG

Lot: A03700

Each vial contains 400ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200µl of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8°C. It is recommended that the reconstituted material be aliquotted and stored at -20°C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

Cross Reaction:

Peptide:	%:
CNP (1-29) (porcine)	100
CNP (1-53) (porcine)	0
CNP (32-53) (porcine)	0
ANF (1-28) (human)	0
BNP-32 (porcine)	0
Endothelin-1 (human, porcine, dog, rat, mouse, bovine)	0
ANF (rat)	0