



Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

T-4245

Rabbit anti-VIP (Vasoactive Intestinal Peptide) (human, mouse, rat) IgG

Lot: A13614

Each vial contains 400ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200µl of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8°C. It is recommended that the reconstituted material be aliquotted and stored at -20°C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

Cross Reaction:

Peptide:	%:
VIP (human, mouse, rat)	100
VIP (guinea pig)	100
VIP (10-28) (human, bovine, porcine, rat)	4.4
VIP (Chicken)	<1
PACAP-27 (human, mouse, ovine, porcine, rat)	<0.02
VIP (1-12) (human, porcine, rat)	0
PHM-27 (human)	0
Substance P	0
Endothelin-1 (human, porcine, dog, rat, mouse, bovine)	0
Secretin (porcine)	0
Glucagon (human)	0
Galanin (rat)	0
Somatostatin-14	0
PACAP-38 amide (human, mouse, ovine, porcine, rat)	0
PACAP-38 (16-38) amide (human, chicken, mouse, ovine, porcine, rat)	0