



Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

T-4271

Rabbit anti-Dynorphin A (1-8) IgG

Lot: A09288

Each vial contains 400ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200 μ l of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8°C. It is recommended that the reconstituted material be aliquotted and stored at -20°C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

Cross Reaction:

Peptide:	%:
Dynorphin A (1-8)	100
Pro-dynorphin (209-240) (porcine)	<0.01
Dynorphin A	<0.01
Dynorphin A (1-13)	<0.01
Dynorphin B	0
β -Endorphin (porcine)	0
Met-Enkephalin	0
α -Neoendorphin (porcine)	0
Leu-Enkephalin-Arg	0