



Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

T-4297

Rabbit anti-Met-Enkephalin-Arg-Gly-Leu IgG

Lot: 990027

Each vial contains 400ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200 μ l of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8°C. It is recommended that the reconstituted material be aliquotted and stored at -20°C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

Cross Reaction:

| Peptide: | %: |
|---------------------------------|-----------|
| Met-Enkephalin-Arg-Gly-Leu | 100 |
| Met-Enkephalin-Arg-Phe | 0 |
| Met-Enkephalin | 0 |
| Leu-Enkephalin | 0 |
| Dynorphin A | 0 |
| β -Endorphin (human) | 0 |
| ACTH (1-39) (human) | 0 |
| α -MSH | 0 |
| Peptide F (bovine) | 0 |
| [Arg ⁸]-Vasopressin | 0 |