



Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

T-4301

Rabbit anti-Met-Enkephalin-Arg-Phe IgG

Lot: 032221

Each vial contains 400ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200µl of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8°C. It is recommended that the reconstituted material be aliquotted and stored at -20°C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

Cross Reaction:

Peptide:	%:
Met-Enkephalin-Arg-Phe	100
Met-Enkephalin-Arg-Gly-Leu	<0.01
Leu-Enkephalin	0
Met-Enkephalin	0
β-Endorphin (human)	0
BAM 22 (bovine)	0
Peptide F (bovine)	0
BAM-3200	0
Dynorphin A	0
ACTH (1-39) (human)	0
Endothelin-1 (human, bovine, porcine, dog, rat, mouse)	0