



Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

T-4390

Rabbit anti-rec IFG (Insulin-like Growth Factor)-1 (human) IgG

Lot: 015769

Each vial contains 400ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200 μ l of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8°C. It is recommended that the reconstituted material be aliquotted and stored at -20°C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

Cross Reaction:

Peptide:	%:
rec IGF-I (human)	100
IGF-I (24-41)	0
rec IGF-II (human)	<0.02
IGF-II (33-40)	0
Insulin (human)	0
Proinsulin	0
Amylin (free acid) (human)	0
Somatostatin-14	0
PACAP-38 (human, mouse, ovine, porcine, rat)	0
Growth Hormone (human)	0
EGF (human)	0