



Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

T-4421

Rabbit anti-Leumorphin (human) IgG

Lot:

Each vial contains 400ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200µl of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8°C. It is recommended that the reconstituted material be aliquotted and stored at -20°C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

Cross Reaction:

Peptide:	%:
Leumorphin (human)	100
Prodynorphin (228-256)	100
Big Dynorphin	28
β-Neoendothelin	0
Dynorphin A (1-8)	0
Dynorphin A (1-17)	0
α-Neoendorphin	0
Leu-Enkephalin	0
Met-Enkephalin	0
Beta-Endorphin (human)	0