



Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

T-4429

Rabbit anti-Mating Factor- α IgG

Lot: A13948

Each vial contains 400ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200 μ l of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8°C. It is recommended that the reconstituted material be aliquotted and stored at -20°C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

Cross Reaction:

| Peptide: | %: |
|-----------------------------------|-----------|
| Mating Factor α | 100 |
| Met-Enkephalin | <0.001 |
| NPY (porcine) | <0.001 |
| Substance P | 0 |
| VIP (human, bovine, porcine, rat) | 0 |
| Somatostatin-14 | 0 |
| CCK Octapeptide (desulfated) | 0 |