



Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

T-4550

Rabbit anti- TGF (Transforming Growth Factor)- α (1-50) (rat) IgG

Lot: A09150

Each vial contains 400ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200 μ l of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8 $^{\circ}$ C. It is recommended that the reconstituted material be aliquotted and stored at -20 $^{\circ}$ C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

Cross Reaction:

Peptide:	%:
TGF- α (1-50) (rat)	100
TGF- α (human)	100
TGF- α (34-43) (rat)	0
rec EGF	0
EGF Receptor Peptide (985-996) (human)	0
GRF (Rat)	0
α -MSH	0
Substance P	0