



Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

T-4554

Rabbit anti- Thr-Ala-Pro-Arg-Atrial Natriuretic Factor (1-28) (Urodilatin) (human, bovine, porcine) - IgG

Lot:

Each vial contains 400ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200µl of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8°C. It is recommended that the reconstituted material be aliquotted and stored at -20°C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

Cross Reaction:

Peptide:	%:
Thr-Ala-Pro-Arg-ANF (1-28) (human)	100
ANF (1-28) (human)	100
β-ANF (dimmer, antiparallel) (human)	100
ANF (rat)	100
BNP-32 (human)	0
CNP (32-53) (porcine, rat, human)	0
Endothelin-1 (human, porcine, canine, rat, mouse, bovine)	0