



Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

T-4737

Rabbit anti- Urotensin II (Human) IgG

Lot: 016024

Each vial contains 400ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200µl of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8°C. It is recommended that the reconstituted material be aliquotted and stored at -20°C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

Cross Reaction:

Peptide:	%:
Urotensin II (human)	100
Urotensin II (goby)	100
Urotensin II (mouse)	10
Urotensin II (rat)	0
Urotensin I	0
CRF (human, rat)	0
Urocortin (human)	0
Urocortin (rat)	0
Somatostatin-14	0
[Arg8]- Vasopressin	0
Angiotensin II	0
Calcitonin (salmon I)	0
NPY (human, rat)	0