



## Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

**T-4775**

### **Rabbit anti- Kisspeptin-13 (4-13) (human) IgG**

**Lot: 019674**

Each vial contains 400ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200µl of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8°C. It is recommended that the reconstituted material be aliquotted and stored at -20°C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

#### **Cross Reaction:**

<b>Peptide:</b>	<b>%:</b>
Kisspeptin-13 (4-13) (human)	100
Malignant Melanoma Metastasis-suppressor	
Kiss-1 (94-121) (human)	100
Pyr-Ser-Leu-Arg-Trp-NH <sub>2</sub>	0
Antho-Rwamide II	0