



## Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

### T-4797

## Rabbit anti- Urotensin II-Related Peptide (human, mouse, rat) - IgG

**Lot: 022174**

Each vial contains 400 ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200µl of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8°C. It is recommended that the reconstituted material be aliquotted and stored at -20°C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

#### Cross Reaction:

Peptide	%
Urotensin II-Related Peptide (human, mouse, rat)	100
Urotensin II (goby)	11.7
Urotensin II (human)	9
Urotensin II (mouse)	5.8
(Pyr <sup>110</sup> )-Prepro-Urotensin II (110-123)	3.6
(Pen <sup>5</sup> )-Urotensin II (4-11) (human)	0
Urotensin I	0

Cross-Reactivity data is taken from Radioimmunoassay applications