



Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

Cat. No. T-4803

Rabbit anti-Intermedin (rat) IgG

Lot:

Specifications

Each vial contains 400ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200µl of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8°C. It is recommended that the reconstituted material be aliquotted and stored at -20°C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

Cross Reaction:

Peptide:	%:
Intermedin (rat)	100
α-CGRP (rat)	5
α-CGRP (human)	3.3
Adrenomedullin (human)	2.6
Adrenomedullin (rat)	4

Cross-Reactivity data is taken from Radioimmunoassay applications