



Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

T-4814

Rabbit anti-Neuromedin S (rat) IgG

Lot: 021596

Each vial contains 400 ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200µl of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8°C. It is recommended that the reconstituted material be aliquotted and stored at -20°C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

Cross Reaction:

Peptide:	%:
Neuromedin S (rat)	100
Neuromedin S (human)	15
Neuromedin U-25 (human)	0
Neuromedin U-23 (rat)	0
PrePro-Neuromedin U (104-136) (human)	0
PrePro-Neuromedin S (70-103) (human)	0
Neuropeptide Y (human, rat)	0
Ghrelin (rat)	0