



Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

T-4816

Rabbit anti-Osteopontin (human) IgG

Lot: 021804

Each vial contains 400 ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200 μ l of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8 $^{\circ}$ C. It is recommended that the reconstituted material be aliquotted and stored at -20 $^{\circ}$ C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

Cross Reactivity:

Peptide:	%:
Osteopontin 170-187 (human)	100
recombinant human Osteopontin (OPN)	25
Osteopontin 164-179 (human)	0
Osteopontin 59-74 (human)	0
Osteopontin (mouse)	0