



## Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

### T-4821

## Rabbit anti-Obestatin (rat) IgG

**Lot: 022141**

Each vial contains 400 ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200 $\mu$ l of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8 $^{\circ}$ C. It is recommended that the reconstituted material be aliquotted and stored at -20 $^{\circ}$ C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

#### Cross Reactivity:

Peptide:	%:
<b>Obestatin (rat)</b>	<b>100</b>
Obestatin (human)	0.2
Ghrelin (rat)	0
Neuropeptide Y (human, rat)	0
ACTH (rat)	0
BNP-32 (rat)	0
$\beta$ -Endorphin (rat)	0