



## Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

### T-4835

## Rabbit anti-Hemokinin (mouse, rat) IgG

**Lot: 022578**

Each vial contains 400 ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200µl of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8°C. It is recommended that the reconstituted material be aliquotted and stored at -20°C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

#### Cross Reactivity:

Peptide:	%:
<b>Hemokinin 1 (mouse, rat)</b>	<b>100</b>
Hemokinin 1 (human)	0.01
Substance P	0.01
Neurokinin A	0
Neurokinin B	0
Neurotensin	0
Oxytocin	0
CRF (human, rat)	0