



Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

T-4844

Rabbit anti- Exendin-4 IgG

Lot: A12888

Each vial contains 400 ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200µl of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8°C. It is recommended that the reconstituted material be aliquotted and stored at -20°C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

Cross Reaction:

PEPTIDE:	%:
Exendin-4	100
Exendin-3	100
Exendin-4 (3-39)	100
Exendin-4 (9-39)	100
GLP-1 (7-36) amide (human, rat, mouse, bovine, guinea pig)	0
GLP-1 (7-37) (human, rat, mouse, bovine, guinea pig))	0
GLP-1 (1-37)(human, rat, mouse, bovine, guinea pig))	0
GLP-2 (1-34)(human)	0
Glucagon (1-29) (human)	0