



Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

T-4846

Rabbit anti- Inactive Ghrelin (human) IgG

Lot: A00474

Each vial contains 400 ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200µl of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8°C. It is recommended that the reconstituted material be aliquotted and stored at -20°C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

Cross Reaction:

PEPTIDE:	%:
(Des-octanoyl)-Ghrelin (human)	100
Ghrelin (human)	0
Motilin (canine)	0
GRF (1-40) (human)	0
GRF (1-43) (rat)	0
Orexin A (human, rat, mouse, bovine)	0
Orexin B (human)	0
Orexin B (rat, mouse)	0
Secretin (human)	0
VIP (human, bovine, porcine, rat)	0
Galanin (human)	0
Galanin (rat)	0