



Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

T-4847

Rabbit anti- Thymosin beta4 (human, bovine, horse, rat) - IgG

Lot: 121207

Each vial contains 400 ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200 μ l of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8 $^{\circ}$ C. It is recommended that the reconstituted material be aliquotted and stored at -20 $^{\circ}$ C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

Cross Reaction:

PEPTIDE	%
Thymosin β -4 (human, bovine, horse, rat)	100
Thymosin β -10	0
Thymosin α -1	0