



Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

Cat. No. T-4852 **Rabbit anti-Neuromedin B IgG**

Lot: A02615

Specifications

Each vial contains 400ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200µl of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8°C. It is recommended that the reconstituted material be aliquotted and stored at -20°C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

CROSS REACTION

PEPTIDE:	%:
Neuromedin B	100%
Neuromedin U-23 (rat)	0%
Neurokinin B	0%
Neurokinin A	0%
Neuromedin U-25 (human)	0%
Neurotensin	0%

* Crossreactivity taken from Enzyme immunoassay application.