



Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

Cat. No. T-4856 **Rabbit anti-Thymosin alpha 1 IgG**

Lot: A10667

Specifications

Each vial contains 400ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200 μ l of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8 $^{\circ}$ C. It is recommended that the reconstituted material be aliquotted and stored at -20 $^{\circ}$ C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

CROSS REACTION

| PEPTIDE: | %: |
|------------------|-----------|
| Thymosin alpha 1 | 100% |
| Thymosin beta 4 | 0% |
| Thymosin beta 10 | 0% |

* Crossreactivity taken from Enzyme immunoassay application.