



## Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

### **Cat. No. T-4858** **Rabbit anti-Spexin IgG**

Lot: A01986

#### **Specifications**

Each vial contains 400ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200µl of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8°C. It is recommended that the reconstituted material be aliquotted and stored at -20°C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

#### **CROSS REACTION**

<b>PEPTIDE:</b>	<b>%:</b>
<b>Spexin</b>	<b>100</b>

\* Crossreactivity taken from Enzyme immunoassay application.