



## Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

### **Cat. No. T-4866** **Rabbit anti-Hepcidin-25 (mouse)**

**Lot# A09277**

#### **Specifications**

Each vial contains 400ug of rabbit IgG in lyophilized form, which was purified from whole serum by protein A column chromatography. Upon reconstitution with 200µl of distilled water, the solution contains 0.1M PBS (pH 7.4) and 0.02% thimerosal.

Lyophilized products are stable for long periods when stored at 2-8°C. It is recommended that the reconstituted material be aliquotted and stored at -20°C (or below) for long-term storage. Repeated freezing and thawing should be avoided.

#### **CROSS REACTION**

<b>PEPTIDE:</b>	<b>%:</b>
<b>Hepcidin-25 (mouse)</b>	<b>100</b>
Hepcidin-25 (rat)	0

\* Crossreactivity taken from Enzyme immunoassay application.