



Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

Cat. No. T-4871

**Rabbit anti-Davunetide [NAP, Activity dependent
Neuroprotective Protein (74-81)] (mouse, rat)**

GENERAL SPECIFICATIONS

STORAGE:

This polyclonal antibody is supplied as a lyophilized powder, suitable for use in EIA/ELISA. For long-term storage, store at -20°C and can be stable for up to one year. Avoid repeated freeze/thaw cycles.

RECONSTITUTION INFORMATION:

Each vial contains enough antiserum for one 96-well plate. The powder should be rehydrated with 5ml of EIA buffer (Catalog # Y-1055).

CROSS REACTION

PEPTIDE:	%:
Davunetide	100