



**Peninsula Laboratories International, Inc.**

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

**Cat. No. Y-1055**  
**EIA Buffer Concentrate**

**Reconstitution Information:**

Each bottle of the EIA Buffer Concentrate contains 50 ml of 20x concentrated EIA buffer. *Dilute the EIA Buffer Concentrate to 1,000 ml (1 liter) with sterile deionized or USP pure water (18 MOhm) and mix well.*

**Composition**

The Buffer Concentrate contains Trizma base and Trizma hydrochloride (pH 7.4), NaCl, Tween 20, 0.2% Thimerosal and 1% BSA.

**Recommended Use**

The EIA Buffer is used for several steps in our EIA (Enzyme Immunoassay) protocols.

**Storage**

Please refrigerate (0 - 4° C).

**Further Information**

Please contact our technical service at: [immunology@bachem.com](mailto:immunology@bachem.com)