



Peninsula Laboratories International, Inc.

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

Cat. No. Y-1060
Standard Diluent-H

Reconstitution Information

Ready to use (8 ml are provided). No dilution required.

Composition

Difibrinated/delipidized charcoal-stripped human serum (Note: no thimerosal or any other preservatives added).

Recommended Use

The Standard Diluent-H is used for standard dilutions as well as for dilutions of human serum or plasma samples.

Storage

Store at -20°C.

Further Information

Please contact our technical service at: immunology@bachem.com.