



**Peninsula Laboratories International, Inc.**

305 Old County Road, San Carlos, CA 94070

Tel: (650) 801-6090

Fax: (650) 595-4071

**Cat. No. Y-1065**  
**Standard Diluent-R**

**Reconstitution Information:**

Ready to use (8 ml are provided). No dilution required.

**Composition**

Charcoal-stripped trace-hemolyzed rat serum. (Note: no thimerosal or any other preservatives added).

**Recommended Use**

The Standard Diluent-R is used for standard dilutions as well as for dilutions of rat serum or plasma samples.

**Storage**

Store at -20°C

**Further Information**

Please contact our technical service at: [immunology@bachem.com](mailto:immunology@bachem.com).